


Pigtail Catheter Sets Hydrophilic




One Step Pigtail with Safety String Lock				
Code	French Size	Length (In Open Position)	Accepts Guide Wire Diameter	Flexible Stiffening Cannula
NC-06191	6F	20cm	0.035"	-
NC-06241	6F	25cm	0.035"	-
NC-06291	6F	30cm	0.035"	-
NC-06391	6F	40cm	0.035"	-
NC-08191	8F	20cm	0.038"	-
NC-08241	8F	25cm	0.038"	-
NC-08291	8F	30cm	0.038"	-
NC-08391	8F	40cm	0.038"	-
NC-10191	10F	20cm	0.038"	-
NC-10241	10F	25cm	0.038"	-
NC-10291	10F	30cm	0.038"	-
NC-10391	10F	40cm	0.038"	-
NC-12191	12F	20cm	0.038"	-
NC-12241	12F	25cm	0.038"	-
NC-12291	12F	30cm	0.038"	-
NC-12391	12F	40cm	0.038"	-
NC-14191	14F	20cm	0.038"	-
NC-14241	14F	25cm	0.038"	-
NC-14291	14F	30cm	0.038"	-
NC-14391	14F	40cm	0.038"	-
NC-08192	8F	20cm	0.038"	+
NC-08242	8F	25cm	0.038"	+
NC-08292	8F	30cm	0.038"	+
NC-08392	8F	40cm	0.038"	+
NC-10192	10F	20cm	0.038"	+
NC-10242	10F	25cm	0.038"	+
NC-10292	10F	30cm	0.038"	+
NC-10392	10F	40cm	0.038"	+
NC-12192	12F	20cm	0.038"	+
NC-12242	12F	25cm	0.038"	+
NC-12292	12F	30cm	0.038"	+
NC-12392	12F	40cm	0.038"	+
NC-14192	14F	20cm	0.038"	+
NC-14242	14F	25cm	0.038"	+
NC-14292	14F	30cm	0.038"	+
NC-14392	14F	40cm	0.038"	+

- Set Contents:**
- Drainage Catheter
 - Safety String Lock
 - Metal Stiffening Cannula
 - Trocar Needle
 - Split Straightener
 - *Flexible Stiffening Cannula (only for 8-14Fr catheter)

- Customized sets available with the following optional contents:
- Connecting tube with 3 way stopcock
 - Scalpel
 - Drainage bag
 - Securement

Qty/Box: 5

Pigtail Catheter Sets Hydrophilic



Closed Loop Pigtail with Safety String Lock (Seldinger Type)

Code	French Size	Length (In Open Position)	Accepts Guide Wire Diameter
NC-06194	6F	20cm	0.035"
NC-06244	6F	25cm	0.035"
NC-06294	6F	30cm	0.035"
NC-06344	6F	35cm	0.035"
NC-06394	6F	40cm	0.035"
NC-08194	8F	20cm	0.038"
NC-08244	8F	25cm	0.038"
NC-08294	8F	30cm	0.038"
NC-08344	8F	35cm	0.038"
NC-08394	8F	40cm	0.038"
NC-10194	10F	20cm	0.038"
NC-10244	10F	25cm	0.038"
NC-10294	10F	30cm	0.038"
NC-10344	10F	35cm	0.038"
NC-10394	10F	40cm	0.038"
NC-12194	12F	20cm	0.038"
NC-12244	12F	25cm	0.038"
NC-12294	12F	30cm	0.038"
NC-12344	12F	35cm	0.038"
NC-12394	12F	40cm	0.038"
NC-14194	14F	20cm	0.038"
NC-14244	14F	25cm	0.038"
NC-14294	14F	30cm	0.038"
NC-14344	14F	35cm	0.038"
NC-14394	14F	40cm	0.038"

Qty/Box: 5


Set Contents:

- Closed Loop Pigtail Drainage Catheter
- Safety String Lock
- Metal Stiffening Cannula
- Flexible Stiffening Cannula
- *Flexible Stiffening Cannula (only for 8-14Fr catheter)

Customized sets available with the following optional contents:

- Connecting tube with 3 way stopcock
- Scalpel
- Drainage bag
- Securement

Pigtail Catheter Sets Hydrophilic



Closed Loop Pigtail (One Step Type)

Code	French Size	Length (In Open Position)	Accepts Guide Wire Diameter	Safety String Lock
NC-06195	6F	20cm	0.035"	+
NC-06245	6F	25cm	0.035"	+
NC-06295	6F	30cm	0.035"	+
NC-06395	6F	40cm	0.035"	+
NC-08195	8F	20cm	0.038"	+
NC-08245	8F	25cm	0.038"	+
NC-08295	8F	30cm	0.038"	+
NC-08395	8F	40cm	0.038"	+
NC-10195	10F	20cm	0.038"	+
NC-10245	10F	25cm	0.038"	+
NC-10295	10F	30cm	0.038"	+
NC-10395	10F	40cm	0.038"	+
NC-12195	12F	20cm	0.038"	+
NC-12245	12F	25cm	0.038"	+
NC-12295	12F	30cm	0.038"	+
NC-12395	12F	40cm	0.038"	+
NC-14195	14F	20cm	0.038"	+
NC-14245	14F	25cm	0.038"	+
NC-14295	14F	30cm	0.038"	+
NC-14395	14F	40cm	0.038"	+
NC-06196	6F	20cm	0.035"	-
NC-06246	6F	25cm	0.035"	-
NC-06296	6F	30cm	0.035"	-
NC-06396	6F	40cm	0.035"	-
NC-08196	8F	20cm	0.038"	-
NC-08246	8F	25cm	0.038"	-
NC-08296	8F	30cm	0.038"	-
NC-08396	8F	40cm	0.038"	-
NC-10196	10F	20cm	0.038"	-
NC-10246	10F	25cm	0.038"	-
NC-10296	10F	30cm	0.038"	-
NC-10396	10F	40cm	0.038"	-
NC-12196	12F	20cm	0.038"	-
NC-12246	12F	25cm	0.038"	-
NC-12296	12F	30cm	0.038"	-
NC-12396	12F	40cm	0.038"	-
NC-14196	14F	20cm	0.038"	-
NC-14246	14F	25cm	0.038"	-
NC-14296	14F	30cm	0.038"	-
NC-14396	14F	40cm	0.038"	-

Set Contents:


- Closed Loop Pigtail Drainage Catheter
- Metal Stiffening Cannula
- Trocar Needle
- Split Straightener

Customized sets available with the following optional contents:

- Connecting tube with 3 way stopcock
- Scalpel
- Drainage bag
- Securement

Qty/Box: 5

Biliary Catheter Sets Hydrophilic



Biliary Catheter with Radiopaque Band and Safety String Lock

Code	French Size	Length (In Open Position)	Accepts Guide Wire Diameter	Flexible Stiffening Cannula
NC-0619B	6F	20cm	0.035"	-
NC-0624B	6F	25cm	0.035"	-
NC-0629B	6F	30cm	0.035"	-
NC-0639B	6F	40cm	0.035"	-
NC-0644B	6F	45cm	0.035"	-
NC-0819B	8F	20cm	0.038"	+
NC-0824B	8F	25cm	0.038"	+
NC-0829B	8F	30cm	0.038"	+
NC-0839B	8F	40cm	0.038"	+
NC-0844B	8F	45cm	0.038"	+
NC-1019B	10F	20cm	0.038"	+
NC-1024B	10F	25cm	0.038"	+
NC-1029B	10F	30cm	0.038"	+
NC-1039B	10F	40cm	0.038"	+
NC-1044B	10F	45cm	0.038"	+
NC-1219B	12F	20cm	0.038"	+
NC-1224B	12F	25cm	0.038"	+
NC-1229B	12F	30cm	0.038"	+
NC-1239B	12F	40cm	0.038"	+
NC-1244B	12F	45cm	0.038"	+
NC-1419B	14F	20cm	0.038"	+
NC-1424B	14F	25cm	0.038"	+
NC-1429B	14F	30cm	0.038"	+
NC-1439B	14F	40cm	0.038"	+
NC-1444B	14F	45cm	0.038"	+

Qty/Box: 5


Set Contents:

- Biliary Drainage Catheter with Radiopaque Band
 - Safety String Lock
 - Metal Stiffening Cannula
 - Flexible Stiffening Cannula
 - Split Straightener
 - MLL Injection Site
- *Flexible Stiffening Cannula (only for 8-14Fr catheter)

Customized sets available with the following optional contents:

- Connecting tube with 3 way stopcock
- Scalpel
- Drainage bag
- Securement


Drainage Catheter Insertion Set



Code	Description	Qty/Box
HD-N6410	Drainage Catheter Insertion Set	5

Set Contents:

- 21Gx15cm Access Needle • OD: 0.5 mm Stylet • 0.018"x 65cm Special Kink Resistant Guidewire
- 22 Gauge Stiffening Cannula • 4Fr x 20cm Dilator • 6Fx18cm Sheath with Hydrophilic Coating



Code	Description	Qty/Box
HD-N6411	Drainage Catheter Insertion Set with 0.035"x150cm Amplatz Guidewire	5
HD-N6412	Drainage Catheter Insertion Set with 0.035"x80cm Amplatz Guidewire	5

Set Contents:

- 21Gx15cm Access Needle • OD: 0.5 mm Stylet • 0.018"x 65cm Special Kink Resistant Guidewire
- 22 Gauge Stiffening Cannula • 4Fr x 20cm Dilator • 6Fx18cm Sheath with Hydrophilic Coating
- 0.035"x 80/150cm Amplatz Guidewire

10mm from tip Radiopaque Band on Sheath

Drainage Catheter Connectors

Code	Description	Stopcock	Qty/Box
DQ-0305	5cm	-	50
DQ-0325	25cm	-	50
DQ-0430	40cm	+	100

