


Glaucus RX PTA Balloon Catheter




*Biometrix RX PTA balloon catheter provides
the best solution for peripheral angioplasty
procedures*

Glaucus RX PTA Balloon Catheter


Biometrix RX PTA balloon catheter is semi-compliant. It is used for dilatation of clogged peripheral blood vessels in angioplasty procedures. The catheter has excellent trackability features and a flexible tapered tip. With optimal pushability and an innovative hydrophilic coating (FC) it is fit for the toughest occlusion crossings.





		Pressure (bars)							
		2	4	6*	8	10	12	14**	16**
Balloon diameter (mm)	2	1,88	1,94	2,00	2,03	2,06	2,09	2,12	2,15
	2,25	2,12	2,18	2,25	2,28	2,32	2,35	2,39	2,42
	2,5	2,36	2,43	2,50	2,54	2,58	2,61	2,65	2,69
	2,75	2,59	2,67	2,75	2,79	2,83	2,87	2,92	2,96
	3	2,83	2,91	3,00	3,05	3,09	3,14	3,18	3,23
	3,5	3,30	3,40	3,50	3,55	3,61	3,66	3,71	3,76
	4	3,77	3,88	4,00	4,06	4,12	4,18	4,24	4,30
	4,5	4,24	4,37	4,50	4,57	4,64	4,70	4,77	4,84
	5	4,80	4,90	5,00	5,08	5,15	5,23	5,30	5,38
	6	5,76	5,88	6,00	6,09	6,18	6,27	6,36	6,45
7	6,72	6,86	7,00	7,11	7,21	7,32	7,42	7,53	

*Nominal pressure **Rated Burst Pressure



ORDERING INFORMATION

Code	Balloon Diameter	Balloon Length
BG-G200008-GLAU	2mm	8mm
BG-G200010-GLAU	2mm	10mm
BG-G200020-GLAU	2mm	20mm
BG-G200030-GLAU	2mm	30mm
BG-G200040-GLAU	2mm	40mm
BG-G200060-GLAU	2mm	60mm
BG-G200080-GLAU	2mm	80mm
BG-G200100-GLAU	2mm	100mm
BG-G200120-GLAU	2mm	120mm
BG-G200150-GLAU	2mm	150mm
BG-G200200-GLAU	2mm	200mm
BG-G225008-GLAU	2,25mm	8mm
BG-G225010-GLAU	2,25mm	10mm
BG-G225020-GLAU	2,25mm	20mm
BG-G225030-GLAU	2,25mm	30mm
BG-G225040-GLAU	2,25mm	40mm
BG-G225060-GLAU	2,25mm	60mm
BG-G225080-GLAU	2,25mm	80mm
BG-G225100-GLAU	2,25mm	100mm
BG-G225120-GLAU	2,25mm	120mm
BG-G225150-GLAU	2,25mm	150mm
BG-G225200-GLAU	2,25mm	200mm
BG-G250008-GLAU	2,5mm	8mm
BG-G250010-GLAU	2,5mm	10mm
BG-G250020-GLAU	2,5mm	20mm
BG-G250030-GLAU	2,5mm	30mm
BG-G250040-GLAU	2,5mm	40mm
BG-G250060-GLAU	2,5mm	60mm
BG-G250080-GLAU	2,5mm	80mm
BG-G250100-GLAU	2,5mm	100mm
BG-G250120-GLAU	2,5mm	120mm
BG-G250150-GLAU	2,5mm	150mm
BG-G250200-GLAU	2,5mm	200mm
BG-G275008-GLAU	2,75mm	8mm
BG-G275010-GLAU	2,75mm	10mm
BG-G275020-GLAU	2,75mm	20mm
BG-G275030-GLAU	2,75mm	30mm
BG-G275040-GLAU	2,75mm	40mm
BG-G275060-GLAU	2,75mm	60mm
BG-G275080-GLAU	2,75mm	80mm
BG-G275100-GLAU	2,75mm	100mm
BG-G275120-GLAU	2,75mm	120mm
BG-G275150-GLAU	2,75mm	150mm

Code	Balloon Diameter	Balloon Length
BG-G300008-GLAU	3mm	8mm
BG-G300010-GLAU	3mm	10mm
BG-G300020-GLAU	3mm	20mm
BG-G300030-GLAU	3mm	30mm
BG-G300040-GLAU	3mm	40mm
BG-G300060-GLAU	3mm	60mm
BG-G300080-GLAU	3mm	80mm
BG-G300100-GLAU	3mm	100mm
BG-G300120-GLAU	3mm	120mm
BG-G300150-GLAU	3mm	150mm
BG-G350008-GLAU	3,50mm	8mm
BG-G350010-GLAU	3,50mm	10mm
BG-G350020-GLAU	3,50mm	20mm
BG-G350030-GLAU	3,50mm	30mm
BG-G350040-GLAU	3,50mm	40mm
BG-G350060-GLAU	3,50mm	60mm
BG-G350080-GLAU	3,50mm	80mm
BG-G350100-GLAU	3,50mm	100mm
BG-G350120-GLAU	3,50mm	120mm
BG-G350150-GLAU	3,50mm	150mm
BG-G400008-GLAU	4mm	8mm
BG-G400010-GLAU	4mm	10mm
BG-G400020-GLAU	4mm	20mm
BG-G400030-GLAU	4mm	30mm
BG-G400040-GLAU	4mm	40mm
BG-G400060-GLAU	4mm	60mm
BG-G400080-GLAU	4mm	80mm
BG-G400100-GLAU	4mm	100mm
BG-G400120-GLAU	4mm	120mm
BG-G450010-GLAU	4,50mm	10mm
BG-G450020-GLAU	4,50mm	20mm
BG-G450030-GLAU	4,50mm	30mm
BG-G450040-GLAU	4,50mm	40mm
BG-G450060-GLAU	4,50mm	60mm
BG-G450080-GLAU	4,50mm	80mm
BG-G600010-GLAU	6mm	10mm
BG-G600020-GLAU	6mm	20mm
BG-G600030-GLAU	6mm	30mm
BG-G600040-GLAU	6mm	40mm
BG-G700010-GLAU	7mm	10mm
BG-G700020-GLAU	7mm	20mm
BG-G700030-GLAU	7mm	30mm
BG-G700040-GLAU	7mm	40mm